Lifetime Limited Warranty

If you have a question about your MOBB Healthcare device or this warranty, please contact an authorized MOBB Healthcare dealer.

MOBB Healthcare warrants to the original consumer purchaser for as long as they own their MOBB product that it will be free from defects in material and workmanship under normal use. This warranty is only applicable to products used in residential applications within the United States and Canada and is not transferable to subsequent owners. Should defects in material or workmanship exist during the applicable warranty period in any MOBB Healthcare Product, MOBB Healthcare will elect at its discretion to either repair or replace the defective material. This warranty covers only parts and materials of the products supplied by MOBB Healthcare. Not covered under this warranty are wear and tear items or abuse of the product in a manner that was not intended for the product to be used. This warranty shall become void if the product is in any way modified, improperly adjusted or abused.

MHADJRTS

990000001315

MOBB HEALTH CARE

4” Adjustable Raised Toilet Seat

MOBB Healthcare Ltd.

116A Viceroy Rd. Unit 5-6
Concord, Ontario, Canada
L4K 2M1
www.mobbhhc.com

Learn about the parts

Assembly Instructions

STEP 1. Remove the item from the box. If the 4" Adjustable Raised Toilet Seat has been exposed to cold, allow it to reach room temperature before assembly. This will help prevent potential cracking during the assembly process.

STEP 2. Lift up the toilet seat and lid.

STEP 3. Place the adjustable toilet seat onto the toilet seat frame, as shown in Figure 2 below.

Assembly Instructions

FIGURE 2

STEP 4.

a. Gently squeeze the sides of the adjustable seat to clear the rim of the toilet seat frame.

Avoid over-squeezing to prevent damaging the seat.

FIGURE 3

Assembly Instructions

STEP 4.
b. Ensure the flange is positioned inside the rim of the toilet seat.
c. Release your grip. The flange will press against the inside of the rim, securing the seat in place.

STEP 5. Confirm that the adjustable seat is securely in place on the toilet seat frame.

CARE & CLEANING

- Regular Cleaning is recommended.
- Use a non-abrasive cleaner or detergent. Rinse with warm water and dry with completely.
- Regularly inspect your adjustable toilet seat for any corrosion, breakage, compression or wear.

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